

signo ortho®
tecnologia a
serviço da vida



SISTEMA FECHAMENTO
DE CRÂNIO 1.5MM

SKULL CLOSURE SYSTEM 1.5MM
SISTEMA CIERRE DE CRÁNEO 1.5MM

SISTEMA FECHAMENTO DE CRÂNIO 1.5MM

SKULL CLOSURE SYSTEM 1.5MM | SISTEMA CIERRE DE CRÁNEO 1.5MM

PARAFUSOS SCREW | TORNILLO


PARAFUSOS 1.5 SCREW 1.5 TORNILLO 1.5			PARAFUSOS 1.5 SCREW 1.5 TORNILLO 1.5		
▼ AUTO PERFORANTE	▼ AUTO ROSQUEANTE	EMERGENCIA 1.8 EMERGENCIA 1.8	▼ AUTO PERFORANTE	▼ AUTO ROSQUEANTE	EMERGENCIA 1.8 EMERGENCIA 1.8
4.0	16.192	16.106	6.0	16.193	16.108
5.0	16.115	16.107	7.0		16.131

▼ AUTO PERFORANTE AUTO DRILLING | AUTO DE PERFORACION
▼ AUTO ROSQUEANTE THREADING AUTO | AUTO ROSQUEANTE


PLACAS (escala 1:1) PLATES (SCALE 1:1) | PLACAS (ESCALA 1:1)

PLACAS 1.5 PLATES 1.5 PLACAS 1.5	ULP #0.3mm	LP #0.6mm
------------------------------------	------------	-----------


Y Y Y #0.6mm	Duplo Y Double Y Doble Y #0.6mm	Duplo T Double T Doble T #0.6mm
----------------	--------------------------------------	--------------------------------------



X X X #0.6mm	Triângulo Triangle Triángulo #0.6mm	Matriz Matrix Matrix #0.6mm
----------------	--	----------------------------------




Reta Straight | Recta #0.6mm




PLACAS (escala 1:1) PLATES (SCALE 1:1) | PLACAS (ESCALA 1:1)

Tela Screen | Pantalla 85x50 #0.2mm | 19.599 203x127 #0.6mm | 19.412



Malha Cruz Cross Stitch | Malla Cruzada 50x50 #0.3mm | 19.641



M DULO DE FECHAMENTO DE CRÂNIO 1.5MM SKULL CLOSURE MODULE 1.5MM | MÓDULO DE CIERRE DE CRÁNEO 1.5MM

Brocas com stop, Pontas de chave, Placas e Parafusos para Fechamento de Crânio 1.5mm

Drills with stop, Drivers' tips, Plates and Screws for Skull Lock 1.5mm

Fresas con stop, Puntas de llave, Placas y Tornillos de Bloqueo de Cráneo 1.5mm



SIGNO ORTHO | SIGNO VINCES

0800 600 5999 | +55 41 3032.5999 | signoortho@signoortho.com.br | www.signoortho.com.br
Rod. Raul Azevedo de Macedo, 733 | CEP 83.606-482 | Campo Largo / PR | Brasil